

TECHNICAL DATA SHEET

Medicom SafeGauze HemoStyp Topical Hemostatic Dressing Items # 4900-STAT and 4900-STATC

Intended Use

Used in Health Care facilities for management of topical wounds and temporary control of external surface bleeding such as minor cuts and minor laceration.

Characteristic	Specification				
Description	Soft, glossy, slightly stretchable woven material with cut and/or weaved edges.				
Material Composition	Carboxymethyl cellulose				
Ply	4 ply				
Color	White or light yellow				
Dimension	Size	Width (cm)		Length (cm)	
	¾" X ¾"	1,9		1,9	
Packaging	Dressing/ Blister	Blisters/ Blister Pack	Blister Packs/ Box	Boxes/ Case	Dressings/ Case
	1	10	2	200	4000

Item Code	Description	Color	Specific Attributes
4900-STAT	Medicom SafeGauze HemoStyp Topical Hemostatic Dressing, 4 ply, 1,9 x 1,9 cm (0.75" x 0.75"), 1 Dressing/Blister, 10 Blisters/Blister Pack, 2 Blister Packs/Box, 200 Boxes/Case, 4000 Dressings/Case	White or Light Yellow	Sterile
4900-STATC	Medicom SafeGauze HemoStyp Topical Hemostatic Dressing, 4 ply, 1,9 x 1,9 cm (0.75" x 0.75"), 1 Dressing/Blister, 10 Blisters/Blister Pack, 2 Blister Packs/Box, 200 Boxes/Case, 4000 Dressings/Case	White or Light Yellow	Sterile

Regulatory

Health Canada Class 2 (Item Code 4900-STATC)

FDA Unclassified, Product Code FRO, 510(k) K112800 (Item Code 4900-STAT)

Storage Conditions

Store in a cool dry place

Country of Origin

China